



# PAXIP1 Antibody

<b>Product Code</b>	CSB-PA017496GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q6ZW49
<b>Immunogen</b>	Human PAXIP1
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB
<b>Storage Buffer</b>	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity purified
<b>Isotype</b>	IgG
<b>Alias</b>	PAX interacting (with transcription-activation domain) protein 1;PAXIP1;CAGF28;CAGF29;FLJ41049;PACIP1;PAXIP1L;PTIP;TNRC2 ;
<b>Product Type</b>	Purified Rabbit Anti human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	PAXIP1
<b>Target Details</b>	This gene is a member of the paired box (PAX) gene family and encodes a nuclear protein with six BRCT (breast cancer carboxy-terminal) domains. This protein plays a critical role in maintaining genome stability, condensation of chromatin and progression through mitosis.
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.