



PDXK Antibody

Product Code	CSB-PA017748GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00764
Immunogen	Human PDXK
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	pyridoxal (pyridoxine, vitamin B6) kinase;PDXK;C21orf124;C21orf97;DKFZp566A071;FLJ31940;FLJ37311;MGC15873;MGC31754;MGC52346;PKH;PNK;PRED79 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PDXK
Target Details	This protein phosphorylates vitamin B6, a step required for the conversion of vitamin B6 to pyridoxal-5-phosphate, an important cofactor in intermediary metabolism. The encoded protein is cytoplasmic and probably acts as a homodimer. Alternatively spliced transcript variants have been described, but their biological validity has not been determined.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.