



PFKFB3 Antibody

Product Code	CSB-PA017819GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16875
Immunogen	Human PFKFB3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	6-phosphofructo-2-kinase/fructose-2,6-biphosphatase 3;PFKFB3;FLJ37326;IPFK2;PFK2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PFKFB3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.