



PHGDH Antibody

Product Code	CSB-PA017920GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43175
Immunogen	Human PHGDH
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	phosphoglycerate dehydrogenase;PHGDH;3-PGDH;3PGDH;MGC3017;PDG;PGAD;PGD;PGDH;SERA ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Metabolism
Target Names	PHGDH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.