



PHYHIP Antibody

Product Code	CSB-PA017948GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92561
Immunogen	Human PHYHIP
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	phytanoyl-CoA 2-hydroxylase interacting protein;PHYHIP;DYRK1AP3;KIAA0273;PAHX-AP;PAHXAP1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PHYHIP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.