



PI16 Antibody

Product Code	CSB-PA017951GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UXB8
Immunogen	Human PI16
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	peptidase inhibitor 16;PI16;CRISP9;DKFZp586B1817;MGC45378;MSMBBP;PSPBP ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PI16
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.