



PIP5K1B Antibody

Product Code	CSB-PA018026GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14986
Immunogen	Human PIP5K1B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	phosphatidylinositol-4-phosphate 5-kinase, type I, beta;PIP5K1B;MSS4;STM7 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PIP5K1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.