



PIP5KL1 Antibody

Product Code	CSB-PA018030GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T9C9
Immunogen	Human PIP5KL1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	phosphatidylinositol-4-phosphate 5-kinase-like 1;PIP5KL1;MGC46424;PIPKH;bA203J24.5 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PIP5KL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.