



PLA2G12A Antibody

Product Code	CSB-PA018086GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZM1
Immunogen	Human PLA2G12A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	phospholipase A2, group X11A;PLA2G12A;GX11;PLA2G12;ROSSY ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PLA2G12A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.