



PLEKHA3 Antibody

Product Code	CSB-PA018156GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HB20
Immunogen	Human PLEKHA3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	pleckstrin homology domain containing, family A (phosphoinositide binding specific) member 3;PLEKHA3;FAPP1;FLJ20067 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PLEKHA3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.