



PLEKHA8 Antibody

Product Code	CSB-PA018162GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96JA3
Immunogen	Human PLEKHA8
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	pleckstrin homology domain containing, family A (phosphoinositide binding specific) member 8;PLEKHA8;FAPP2;MGC3358 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PLEKHA8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.