



PMP2 Antibody

Product Code	CSB-PA018240GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02689
Immunogen	Human PMP2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	peripheral myelin protein 2;PMP2;FABP8;M-FABP;MP2;P2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Transport
Target Names	PMP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.