



PNKD Antibody

Product Code	CSB-PA018260GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N490
Immunogen	Human PNKD
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	paroxysmal nonkinesigenic dyskinesia;PNKD;BRP17;DKFZp564N1362;DYT8;FKSG19;FPD1;KIAA1184;KI PP1184;MGC31943;MR-1;MR1;PDC;TAHCCP2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PNKD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.