



PNRC2 Antibody

Product Code	CSB-PA018290GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPJ4
Immunogen	Human PNRC2
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	proline-rich nuclear receptor coactivator 2;PNRC2;FLJ20312;MGC99541 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PNRC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.