



POLG2 Antibody

Product Code	CSB-PA018314GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHN1
Immunogen	Human POLG2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	polymerase (DNA directed), gamma 2, accessory subunit;POLG2;HP55;MTPOLB;PEOA4;POLB;POLG-BETA;POLGB ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	POLG2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.