



POLR1D Antibody

Product Code	CSB-PA018325GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2S0
Immunogen	Human POLR1D
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	polymerase (RNA) I polypeptide D, 16kDa;POLR1D;FLJ20616;MGC9850;POLR1C;RPA16;RPA9;RPAC2;RPO1-3 ;; hRPA19; AC19
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	POLR1D
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.