



POLR2L Antibody

Product Code	CSB-PA018341GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62875
Immunogen	Human POLR2L
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	polymerase (RNA) II (DNA directed) polypeptide L, 7.6kDa;POLR2L;RBP10;RPABC5;RPB10;RPB10beta;RPB7.6;hRPB7.6;hsRPB 10b ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	POLR2L
Target Details	This gene encodes a subunit of RNA polymerase II, the polymerase responsible for synthesizing messenger RNA in eukaryotes. The product of this gene contains four conserved cysteines characteristic of an atypical zinc-binding domain. Like its counterpart in yeast, this subunit may be shared by the other two DNA-directed RNA polymerases.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.