



POSTN Antibody

Product Code	CSB-PA018381GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15063
Immunogen	Human POSTN
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	periostin, osteoblast specific factor;POSTN;MGC119510;MGC119511;OSF-2;OSF2;PDLPOSTN;PN;RP11-412K4.1;periostin ;
Product Type	Purified Rabbit Anti Homo sapiens PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Adhesion
Target Names	POSTN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.