



PPAPDC3 Antibody

Product Code	CSB-PA018420GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NBV4
Immunogen	Human PPAPDC3
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	phosphatidic acid phosphatase type 2 domain containing 3;PPAPDC3;C9orf67;FLJ14662;KIAA0515;MGC12921;NET39 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PLPP7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.