



PPIG Antibody

Product Code	CSB-PA018477GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13427
Immunogen	Human PPIG
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	peptidylprolyl isomerase G (cyclophilin G);PPIG;CARS-Cyp;CYP;MGC133241;SRCyp ;CASP10
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPIG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.