



PPM1K Antibody

Product Code	CSB-PA018497GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N3J5
Immunogen	Human PPM1K
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	protein phosphatase 1K (PP2C domain containing);PPM1K;DKFZp667B084;DKFZp761G058;PP2Cm;PTMP;UG0882E07 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPM1K
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.