



PPM1M Antibody

Product Code	CSB-PA018499GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96MI6
Immunogen	Human PPM1M
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	protein phosphatase 1M (PP2C domain containing); PPM1M; FLJ32332; PP2C-eta; PP2CE; PP2Ceta ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPM1M
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.