



PPP1R13L Antibody

Product Code	CSB-PA018518GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUF5
Immunogen	Human PPP1R13L
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	protein phosphatase 1, regulatory (inhibitor) subunit 13 like;PPP1R13L;IASPP;NKIP1;RAI ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPP1R13L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.