



PPP1R14A Antibody

Product Code	CSB-PA018519GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96A00
Immunogen	Human PPP1R14A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	protein phosphatase 1, regulatory (inhibitor) subunit 14A;PPP1R14A;CPI-17;CPI17;PPP1INL ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPP1R14A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.