



PPP1R9B Antibody

Product Code	CSB-PA018558GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96SB3
Immunogen	Human PPP1R9B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	0.1M NaHCO ₃ ,0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	protein phosphatase 1, regulatory (inhibitor) subunit 9B;PPP1R9B;FLJ30345;PPP1R6;PPP1R9;SPINO;Spn ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPP1R9B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.