



# PRKDC Antibody

<b>Product Code</b>	CSB-PA018714GA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P78527
<b>Immunogen</b>	Human PRKDC
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	ELISA, WB
<b>Storage Buffer</b>	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
<b>Purification Method</b>	Antigen Affinity purified
<b>Isotype</b>	IgG
<b>Alias</b>	protein kinase, DNA-activated, catalytic polypeptide; PRKDC; DNA-PKcs; DNAPK; DNP1; HYRC; HYRC1; XRCC7; p350 ;
<b>Product Type</b>	Purified Rabbit Anti Human PolyClonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Others
<b>Target Names</b>	PRKDC
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.