



PRND Antibody

Product Code	CSB-PA018738GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKY0
Immunogen	Human Prnd
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	prion protein dublet; Prnd; A1450264; Dpl; PrPLP; doppel ;
Product Type	Purified Rabbit Anti Mus musculus PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRND
Target Details	This gene is found on chromosome 20, approximately 20 kbp downstream of the gene encoding cellular prion protein, to which it is biochemically and structurally similar. This protein is a membrane glycosylphosphatidylinositol-anchored glycoprotein that is found predominantly in testis. Mutations in this gene may lead to neurological disorders.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.