



PSG1 Antibody

Product Code	CSB-PA018849GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11464
Immunogen	Human PSG1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	pregnancy specific beta-1-glycoprotein 1; PSG1; B1G1; CD66f; DHFRP2; FLJ90598; FLJ90654; PBG1; PSBG1; PSGGA; PS GIIA; SP1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PSG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.