



PSMA1 Antibody

Product Code	CSB-PA018865GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25786
Immunogen	Human PSMA1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	proteasome (prosome, macropain) subunit, alpha type, 1;PSMA1;HC2;MGC14542;MGC14575;MGC14751;MGC1667;MGC21459;MGC 22853;MGC23915;NU;PROS30 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PSMA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.