



PSME4 Antibody

Product Code	CSB-PA018919GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14997
Immunogen	Human PSME4
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	proteasome (prosome, macropain) activator subunit 4;PSME4;FLJ21864;KIAA0077;MGC138374;MGC142228;PA200 ;
Product Type	Purified Rabbit Anti Human,Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PSME4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.