



PSMG2 Antibody

Product Code	CSB-PA018922GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969U7
Immunogen	Human PSMG2
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	proteasome (prosome, macropain) assembly chaperone 2; PSMG2; CLAST3; HCCA3; HsT1707; MDS003; MGC15092; PAC2; TNFSF5IP1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PSMG2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.