



PTPN23 Antibody

Product Code	CSB-PA019037GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H3S7
Immunogen	Human PTPN23
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	protein tyrosine phosphatase, non-receptor type 23;PTPN23;DKFZp564F0923;HD-PTP;HDPTP;KIAA1471;PTP-TD14 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PTPN23
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.