



PTPRCAP Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA019050GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q14761 |
| Immunogen | Human PTPRCAP |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, WB |
| Storage Buffer | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Alias | protein tyrosine phosphatase, receptor type, C-associated protein; PTPRCAP; CD45-AP; LPAP; MGC138602; MGC138603 ; |
| Product Type | Purified Rabbit Anti Human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | PTPRCAP |
| Target Details | This protein was identified as a transmembrane phosphoprotein specifically associated with tyrosine phosphatase PTPRC/CD45, a key regulator of T- and B-lymphocyte activation. The interaction with PTPRC may be required for the stable expression of this protein. |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |