



PYCR2 Antibody

Product Code	CSB-PA019116GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96C36
Immunogen	Human PYCR2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	0.1M NaHCO ₃ , 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	pyrroline-5-carboxylate reductase family, member 2; PYCR2; FLJ54750; P5CR2 ;
Product Type	Purified Rabbit Anti Human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PYCR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.