



PZP Antibody

Product Code	CSB-PA019131GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20742
Immunogen	Human PZP
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	pregnancy-zone protein;PZP;CPAMD6;MGC133093 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PZP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.