



RAB40B Antibody

Product Code	CSB-PA019204GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12829
Immunogen	Human RAB40B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	RAB40B, member RAS oncogene family;RAB40B;FLJ42385;RAR;SEC4L ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RAB40B
Target Details	This protein has similarity to a yeast protein which suggests a role of the gene product in regulating secretory vesicles.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.