



RAB6A Antibody

Product Code	CSB-PA019216GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20340
Immunogen	Human RAB6A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3.-20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	RAB6A, member RAS oncogene family;RAB6A;RAB6;RAB6A;RAB6B ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RAB6A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.