



# RACGAP1 Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA019250GA01HU   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                        |
| <b>Uniprot No.</b>         | Q9H0H5   |
| <b>Immunogen</b>           | Human RACGAP1  |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human,Mouse,Rat  |
| <b>Tested Applications</b> | ELISA,WB,IHC,IF  |
| <b>Storage Buffer</b>      | PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| <b>Purification Method</b> | Antigen Affinity Purified  |
| <b>Isotype</b>             | IgG  |
| <b>Alias</b>               | Rac GTPase activating protein 1;RACGAP1;HsCYK-4;ID-GAP;MgcRacGAP ;                   |
| <b>Product Type</b>        | Purified Rabbit Anti human PolyClonal Antibody                                       |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Target Names</b>        | RACGAP1  |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.          |