



RANGRF Antibody

Product Code	CSB-PA019319GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HD47
Immunogen	Human RANGRF
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	RAN guanine nucleotide release factor;RANGRF;DKFZp686F02139;HSPC165;HSPC236;MGC110973;MOG1;RANGNRF ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RANGRF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.