



RAP1B Antibody

Product Code	CSB-PA019323GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61224
Immunogen	Human RAP1B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	RAP1B, member of RAS oncogene family;RAP1B;DKFZp586H0723;K-REV;RAL1B ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RAP1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.