



RBPJ Antibody

Product Code	CSB-PA019486GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06330
Immunogen	Human RBPJ
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	recombination signal binding protein for immunoglobulin kappa J region;RBPJ;CBF1;IGKJRB;IGKJRB1;KBF2;MGC61669;RBP-J;RBPJK;RBPSUH;SUH;cs1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RBPJ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.