



RCAN2 Antibody

Product Code	CSB-PA019501GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14206
Immunogen	Human RCAN2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	regulator of calcineurin 2;RCAN2;CSP2;DSCR1L1;MCIP2;RCN2;ZAKI-4;ZAKI4;hRCN2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RCAN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.