



RDH11 Antibody

Product Code	CSB-PA019523GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TC12
Immunogen	Human RDH11
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	retinol dehydrogenase 11 (all-trans/9-cis/11-cis);RDH11;ARSDR1;CGI-82;FLJ32633;HCBP12;MDT1;PSDR1;RALR1;SCALD ;SDR7C1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RDH11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.