



REEP5 Antibody

Product Code	CSB-PA019544GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00765
Immunogen	Human REEP5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	receptor accessory protein 5;REEP5;C5orf18;D5S346;DP1;MGC70440;TB2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	REEP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.