



RHPN1 Antibody

Product Code	CSB-PA019703GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCX5
Immunogen	Human RHPN1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	rhophilin, Rho GTPase binding protein 1; RHPN1; KIAA1929; ODF5; RHOPHILIN; RHPN ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RHPN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.