



RNMT Antibody

Product Code	CSB-PA019902GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43148
Immunogen	Human RNMT
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	RNA (guanine-7-) methyltransferase;RNMT;DKFZp686H1252;KIAA0398;MET;RG7MT1;hCMT1c ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RNMT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.