



RPA1 Antibody

Product Code	CSB-PA020088GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27694
Immunogen	Human RPA1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	replication protein A1, 70kDa;RPA1;HSSB;MST075;REPA1;RF-A;RP-A;RPA70 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RPA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.