



RUVBL1 Antibody

Product Code	CSB-PA020598GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y265
Immunogen	Human RUVBL1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IP
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	RuvB-like 1 (E. coli);RUVBL1;ECP54;INO80H;NMP238;PONTIN;Pontin52;RVB1;TIH1;TIP49;TI P49A ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RUVBL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.