



S100B Antibody

Product Code	CSB-PA020643GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04271
Immunogen	Human S100B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Form	Liquid
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	S100 calcium binding protein B;S100B;NEF;S100;S100beta ;
Immunogen Species	Homo sapiens (Human)
Target Names	S100B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.