



S1PR5 Antibody

Product Code	CSB-PA020654GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H228
Immunogen	Human S1PR5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	sphingosine-1-phosphate receptor 5;S1PR5;EDG8;Edg-8;S1P5;SPPR-1;SPPR-2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	S1PR5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.